

Product datasheet for **HP200501**

MC2 receptor (MC2R) Human qPCR Primer Pair (NM_000529)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4158
Forward Sequence:	GTGGTGCTTACGGTCATCTGGA
Reverse Sequence:	AGGCACAGGATGAAGACCAGCA
ACCN:	<u>NM_000529</u> , <u>NM_000529.1</u> , <u>NM_000529.2</u> , <u>BC094710</u> , <u>BC094710.1</u> , <u>BC069074</u> , <u>BC104169</u> , <u>BC104170</u> , <u>BP196054</u>
UniProt ID:	<u>Q01718</u>
Synonyms:	ACTHR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)