

## Product datasheet for **HP200503**

### Ornithine Carbamoyltransferase (OTC) Human qPCR Primer Pair (NM\_000531)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5009
Forward Sequence:	ATCCTGGCTGATTACCTCACGC
Reverse Sequence:	CCGAATTCGCTGCGCTCATCA
ACCN:	<a href="#">NM_000531</a> , <a href="#">NM_000531.1</a> , <a href="#">NM_000531.2</a> , <a href="#">NM_000531.3</a> , <a href="#">NM_000531.4</a> , <a href="#">NM_000531.5</a> , <a href="#">BC074745</a> , <a href="#">BC107153</a> , <a href="#">BC107154</a> , <a href="#">BC114496</a> , <a href="#">BX118743</a> , <a href="#">BX494966</a> , <a href="#">NM_000531.6</a>
UniProt ID:	<a href="#">P00480</a>
Synonyms:	OCTD; OTCD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)