

## Product datasheet for **HP200526**

### **CORD2 (CRX) Human qPCR Primer Pair (NM\_000554)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	1406
<b>Forward Sequence:</b>	CCAGTGTGGATCTGATGCACCA
<b>Reverse Sequence:</b>	GGTACTGGGTCTTGCAACAG
<b>ACCN:</b>	<u><a href="#">NM_000554</a></u> , <u><a href="#">NM_000554.1</a></u> , <u><a href="#">NM_000554.2</a></u> , <u><a href="#">NM_000554.3</a></u> , <u><a href="#">NM_000554.4</a></u> , <u><a href="#">NM_000554.5</a></u> , <u><a href="#">BC016664</a></u> , <u><a href="#">BC016664.1</a></u> , <u><a href="#">BC009045</a></u> , <u><a href="#">BC053672</a></u> , <u><a href="#">BE782215</a></u> , <u><a href="#">BG396702</a></u> , <u><a href="#">BQ423938</a></u> , <u><a href="#">BT007364</a></u> , <u><a href="#">BU736672</a></u> , <u><a href="#">NM_000554.6</a></u>
<b>UniProt ID:</b>	<u><a href="#">O43186</a></u>
<b>Synonyms:</b>	CORD2; CRD; LCA7; OTX3
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)