

Product datasheet for **HP200531**

CD53 Human qPCR Primer Pair (NM_000560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	963
Forward Sequence:	GCTATGCGAAAGCAAGACTGTGG
Reverse Sequence:	GTCAATCTGGCAGTTCAGGGTC
ACCN:	NM_000560 , NM_000560.1 , NM_000560.2 , NM_000560.3 , BC040693 , BC040693.1 , BC035456 , BG536920 , BT019643 , NM_000560.4
UniProt ID:	P19397
Synonyms:	MOX44; TSPAN25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)