

Product datasheet for **HP200532**

GSTM1 Human qPCR Primer Pair (NM_000561)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2944
Forward Sequence:	TGATGTCCTTGACCTCCACCGT
Reverse Sequence:	GCTGGACTTCATGTAGGCAGAG
ACCN:	<u>NM_000561</u> , <u>NM_000561.1</u> , <u>NM_000561.2</u> , <u>NM_000561.3</u> , <u>BC024005</u> , <u>BC024005.2</u> , <u>BC036805</u> , <u>BC095463</u> , <u>BQ436476</u> , <u>NM_000561.4</u>
UniProt ID:	<u>P09488</u>
Synonyms:	GST1; GSTM1-1; GSTM1a-1a; GSTM1b-1b; GTH4; GTM1; H-B; MU; MU-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)