

Product datasheet for **HP200535**

IL6R Human qPCR Primer Pair (NM_000565)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3570
Forward Sequence:	GACTGTGCACTTGCTGGTGGAT
Reverse Sequence:	ACTTCCTACCAAGAGCACAGC
ACCN:	<u>NM_000565</u> , <u>NM_000565.1</u> , <u>NM_000565.2</u> , <u>NM_000565.3</u> , <u>BC089410</u> , <u>BC035327</u> , <u>BC132684</u> , <u>BC132686</u> , <u>BM911928</u> , <u>NM_000565.4</u>
UniProt ID:	<u>P08887</u>
Synonyms:	CD126; gp80; HIES5; IL-1Ra; IL-6R; IL-6R-1; IL-6RA; IL6Q; IL6QTL; IL6RA; IL6RQ
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)