

Product datasheet for **HP200537**

C Reactive Protein (CRP) Human qPCR Primer Pair (NM_000567)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1401
Forward Sequence:	TCGTGGAGTTCTGGGTAGATGG
Reverse Sequence:	TTCCACCGAAGGAATCCTGCT
ACCN:	NM_000567 , NM_000567.1 , NM_000567.2 , BC125135 , BC020766 , BC070257 , BG533281 , BG618397 , BG618967 , NM_000567.3
UniProt ID:	P02741
Synonyms:	PTX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)