

Product datasheet for **HP200543**

IL1 alpha (IL1A) Human qPCR Primer Pair (NM_000575)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3552
Forward Sequence:	TGTATGTGACTGCCCAAGATGAAG
Reverse Sequence:	AGAGGAGGTTGGTCTCACTACC
ACCN:	<u>BC013142</u> , <u>NM_000575</u> , <u>NM_000575.1</u> , <u>NM_000575.2</u> , <u>NM_000575.3</u> , <u>NM_000575.4</u> , <u>BC013142.2</u> , <u>BC035169</u> , <u>BT007014</u> , <u>NM_000575.5</u>
UniProt ID:	<u>P01583</u>
Synonyms:	IL-1 alpha; IL-1A; IL1; IL1-ALPHA; IL1F1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)