

Product datasheet for **HP200556**

IL4 Human qPCR Primer Pair (NM_000589)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3565
Forward Sequence:	CCGTAACAGACATCTTTGCTGCC
Reverse Sequence:	GAGTGTCTTCTCATGGTGGCT
ACCN:	NM_000589 , NM_000589.1 , NM_000589.2 , NM_000589.3 , BC070123 , BC070123.1 , BC066277 , BC066278 , BC067514 , BC067515 , NM_000589.4
UniProt ID:	P05112
Synonyms:	BCGF-1; BCGF1; BSF-1; BSF1; IL-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)