

## Product datasheet for **HP200567**

### IL6 Human qPCR Primer Pair (NM\_000600)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3569
Forward Sequence:	AGACAGCCACTCACCTCTTCAG
Reverse Sequence:	TTCTGCCAGTGCCTCTTTGCTG
ACCN:	<u><a href="#">BC015511</a></u> , <u><a href="#">NM_000600</a></u> , <u><a href="#">NM_000600.1</a></u> , <u><a href="#">NM_000600.2</a></u> , <u><a href="#">NM_000600.3</a></u> , <u><a href="#">NM_000600.4</a></u> , <u><a href="#">BT019748</a></u> , <u><a href="#">BT019749</a></u> , <u><a href="#">NM_000600.5</a></u>
UniProt ID:	<u><a href="#">P05231</a></u>
Synonyms:	BSF-2; BSF2; CDF; HGF; HSF; IFN-beta-2; IFNB2; IL-6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)