

## Product datasheet for **HP200569**

### PAI1 (SERPINE1) Human qPCR Primer Pair (NM\_000602)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5054
Forward Sequence:	CTCATCAGCCACTGGAAAGGCA
Reverse Sequence:	GACTCGTGAAGTCAGCCTGAAAC
ACCN:	<a href="#">NM_000602</a> , <a href="#">NM_000602.1</a> , <a href="#">NM_000602.2</a> , <a href="#">NM_000602.3</a> , <a href="#">NM_000602.4</a> , <a href="#">BC010860</a> , <a href="#">BC010860.1</a> , <a href="#">BC005927</a> , <a href="#">BC020765</a> , <a href="#">BF351821</a> , <a href="#">NM_000602.5</a>
UniProt ID:	<a href="#">P05121</a>
Synonyms:	PAI; PAI-1; PAI1; PLANH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)