

Product datasheet for **HP200570**

eNOS (NOS3) Human qPCR Primer Pair (NM_000603)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4846
Forward Sequence:	GAAGGCGACAATCCTGTATGGC
Reverse Sequence:	TGTTTCGAGGGACACCACGTCAT
ACCN:	NM_000603 , NM_000603.1 , NM_000603.2 , NM_000603.3 , NM_000603.4 , BC063294 , BC063294.1 , BC069465 , BU731828 , NM_000603.5
UniProt ID:	P29474
Synonyms:	ECNOS; eNOS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)