

Product datasheet for **HP200572**

Interferon alpha 2 (IFNA2) Human qPCR Primer Pair (NM_000605)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3440
Forward Sequence:	TGGGCTGTGATCTGCCTCAAAC
Reverse Sequence:	CAGCCTTTTGGAACTGGTTGCC
ACCN:	<u>NM_000605</u> , <u>NM_000605.1</u> , <u>NM_000605.2</u> , <u>NM_000605.3</u> , <u>BC074936</u> , <u>BC074937</u> , <u>BC104163</u> , <u>BC104164</u> , <u>NM_000605.4</u>
UniProt ID:	<u>P01563</u>
Synonyms:	IFN-alpha-2; IFN-alphaA; IFNA; IFNA2B; IeIF A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)