

## Product datasheet for **HP200574**

### Alpha 1 Acid Glycoprotein (ORM1) Human qPCR Primer Pair (NM\_000607)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5004
Forward Sequence:	CTGACAAGCCAGAGACGACCAA
Reverse Sequence:	TGCTTCTCCAGTGGCTCACACT
ACCN:	<a href="#">NM_000607</a> , <a href="#">NM_000607.1</a> , <a href="#">NM_000607.2</a> , <a href="#">BC026238</a> , <a href="#">BC026238.1</a> , <a href="#">BC007032</a> , <a href="#">BC104818</a> , <a href="#">BC104820</a> , <a href="#">BC143313</a> , <a href="#">BC143314</a> , <a href="#">BI908589</a> , <a href="#">BT019790</a> , <a href="#">NM_000607.4</a>
UniProt ID:	<a href="#">P02763</a>
Synonyms:	AGP-A; AGP1; HEL-S-153w; ORM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)