

Product datasheet for **HP200604**

Fas Ligand (FASLG) Human qPCR Primer Pair (NM_000639)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	356
Forward Sequence:	GGTTCTGGTTGCCTTGGTAGGA
Reverse Sequence:	CTGTGTGCATCTGGCTGGTAGA
ACCN:	<u>NM_000639</u> , <u>NM_000639.1</u> , <u>NM_000639.2</u> , <u>BC017502</u> , <u>BC017502.1</u> , <u>NM_000639.3</u>
UniProt ID:	<u>P48023</u>
Synonyms:	ALPS1B; APT1LG1; APTL; CD95-L; CD95L; CD178; FASL; TNFSF6; TNLG1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)