

Product datasheet for **HP200631**

alpha 2b Adrenergic Receptor (ADRA2B) Human qPCR Primer Pair (NM_000682)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	151
Forward Sequence:	TCTTCGGCAACGCTCTGGTCAT
Reverse Sequence:	CCAGCGAGAAAGGGATGATGAG
ACCN:	NM_000682 , NM_000682.1 , NM_000682.2 , NM_000682.3 , NM_000682.4 , NM_000682.5 , NM_000682.6 , BC133021 , BC136537 , BF475329 , NM_000682.7
UniProt ID:	P18089
Synonyms:	ADRA2L1; ADRA2RL1; ADRARL1; alpha-2BAR; ALPHA2BAR; FAME2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)