

## Product datasheet for **HP200669**

### CACNA1E Human qPCR Primer Pair (NM\_000721)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	777
Forward Sequence:	AGCGTGAGACAGGCAAAGCCAT
Reverse Sequence:	GGATGCACATCTCAAAGTAGCGC
ACCN:	<u>NM_000721</u> , <u>NM_000721.1</u> , <u>NM_000721.2</u> , <u>NM_000721.3</u> , <u>BC037800</u> , <u>BC172376</u> , <u>BE504178</u> , <u>BI914369</u> , <u>BP230397</u> , <u>NM_000721.4</u>
UniProt ID:	<u>Q15878</u>
Synonyms:	BII; CACH6; CACNL1A6; Cav2.3; DEE69; EIEE69; gm139
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)