

Product datasheet for **HP200670**

CACNA2D1 Human qPCR Primer Pair (NM_000722)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	781
Forward Sequence:	GCTATTCACGGATGGAGGAGAAG
Reverse Sequence:	CCATCCACTGAATAGGTCCTCTG
ACCN:	NM_000722 , NM_000722.1 , NM_000722.3 , BC051339 , BC117468 , BC117470 , BG055272 , BI545104 , BP258070 , BU501177 , NM_000722.4
UniProt ID:	P54289
Synonyms:	CACNA2; CACNL2A; CCHL2A; LINC01112; lncRNA-N3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)