

Product datasheet for **HP200684**

CHRM1 Human qPCR Primer Pair (NM_000738)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1128
Forward Sequence:	CTGGCTGGTTTCCTTTGTGCTC
Reverse Sequence:	GGAGAGGAACTGGATGTAGCAC
ACCN:	BC007740 , NM_000738 , NM_000738.1 , NM_000738.2 , BC007740.2 , BC022984 , BC022984.1 , BT007166 , NM_000738.3
UniProt ID:	P11229
Synonyms:	HM1; M1; M1R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)