

Product datasheet for **HP200694**

Nicotinic Acetylcholine Receptor beta 2 (CHRNA2) Human qPCR Primer Pair (NM_000748)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1141
Forward Sequence:	CTCCCTTCCAAACACATCTGGC
Reverse Sequence:	GATGCTGCCATCATAGGAGACC
ACCN:	<u>NM_000748, NM_000748.1, NM_000748.2, BC075040, BC075040.2, BC075041, NM_000748.3</u>
UniProt ID:	<u>P17787</u>
Synonyms:	EFNL3; nAChRB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)