

Product datasheet for **HP200698**

COMT Human qPCR Primer Pair (NM_000754)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1312
Forward Sequence:	CACACTGGACATGGTCTTCCTC
Reverse Sequence:	GCAGATCACGTTGTCAGCCAGT
ACCN:	<u>NM_000754</u> , <u>NM_000754.1</u> , <u>NM_000754.2</u> , <u>NM_000754.3</u> , <u>BC005867</u> , <u>BC005867.2</u> , <u>BC011935</u> , <u>BC011935.2</u> , <u>BC100018</u> , <u>BC000419</u> , <u>BE386545</u> , <u>BE728779</u> , <u>BT007125</u> , <u>BX461394</u> , <u>NM_000754.4</u>
UniProt ID:	<u>P21964</u>
Synonyms:	HEL-S-98n
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)