

Product datasheet for **HP200702**

GM CSF (CSF2) Human qPCR Primer Pair (NM_000758)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1437
Forward Sequence:	GGAGCATGTGAATGCCATCCAG
Reverse Sequence:	CTGGAGGTCAAACATTTCTGAGAT
ACCN:	<u>NM_000758</u> , <u>NM_000758.1</u> , <u>NM_000758.2</u> , <u>NM_000758.3</u> , <u>BC108724</u> , <u>BC108724.1</u> , <u>BC113924</u> , <u>BC113999</u> , <u>NM_000758.4</u>
UniProt ID:	<u>P04141</u>
Synonyms:	CSF; GMCSF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)