

Product datasheet for **HP200710**

Cytochrome P450 2B6 (CYP2B6) Human qPCR Primer Pair (NM_000767)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1555
Forward Sequence:	ACAGTGTGGAGAAGCACCGTGA
Reverse Sequence:	GGTTGAGGTTCTGGTGGCTGAA
ACCN:	<u>NM_000767</u> , <u>NM_000767.2</u> , <u>NM_000767.3</u> , <u>NM_000767.4</u> , <u>BC067430</u> , <u>BC067431</u> , <u>NM_000767.5</u>
UniProt ID:	<u>P20813</u>
Synonyms:	CPB6; CYP2B; CYP2B7; CYP2B7P; CYP11B6; EFVM; IIB1; P450
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)