

Product datasheet for **HP200733**

DIO1 Human qPCR Primer Pair (NM_000792)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1733
Forward Sequence:	CACTGCCTGAGAGGCTCTACAT
Reverse Sequence:	CCAGAACAGCACGAACTTCCTC
ACCN:	NM_000792 , NM_000792.1 , NM_000792.2 , NM_000792.3 , NM_000792.4 , NM_000792.5 , NM_000792.6 , BC017955 , BC107170 , BC107171 , BE046205 , NM_000792.7
UniProt ID:	P49895
Synonyms:	5DI; TXDI1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)