

Product datasheet for **HP200741**

FGF1 Human qPCR Primer Pair (NM_000800)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2246
Forward Sequence:	ATGGCACAGTGGATGGGACAAG
Reverse Sequence:	TAAAAGCCCGTCGGTGCCATG
ACCN:	<u>NM_000800</u> , <u>NM_000800.1</u> , <u>NM_000800.2</u> , <u>NM_000800.3</u> , <u>NM_000800.4</u> , <u>BC032697</u> , <u>BC032697.1</u> , <u>BG706412</u> , <u>NM_000800.5</u>
UniProt ID:	<u>P05230</u>
Synonyms:	AFGF; ECGF; ECGF-beta; ECGFA; ECGFB; FGF-1; FGF-alpha; FGFA; GLIO703; HBGF-1; HBGF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)