

## Product datasheet for **HP200747**

### GABA A Receptor alpha 1 (GABRA1) Human qPCR Primer Pair (NM\_000806)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2554
Forward Sequence:	CACAAGTCTCCTTCTGGCTCAAC
Reverse Sequence:	GGAGTTTCTGGCACTGATGCTC
ACCN:	<a href="#">NM_000806</a> , <a href="#">NM_000806.1</a> , <a href="#">NM_000806.2</a> , <a href="#">NM_000806.3</a> , <a href="#">NM_000806.4</a> , <a href="#">NM_000806.5</a> , <a href="#">BC030696</a> , <a href="#">BC030696.1</a> , <a href="#">BM927274</a>
UniProt ID:	<a href="#">P14867</a>
Synonyms:	DEE19; ECA4; EIEE19; EJM; EJM5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)