

Product datasheet for **HP200755**

GABA A Receptor beta 3 (GABRB3) Human qPCR Primer Pair (NM_000814)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2562
Forward Sequence:	CAGCCAAGGCAAAGAATGACCG
Reverse Sequence:	ATGCCGCCTGAGACCTCATTCA
ACCN:	<u>NM_000814</u> , <u>NM_000814.1</u> , <u>NM_000814.2</u> , <u>NM_000814.3</u> , <u>NM_000814.4</u> , <u>NM_000814.5</u> , <u>BC010641</u> , <u>BC038837</u> , <u>BE220174</u> , <u>BE783170</u> , <u>NM_000814.6</u>
UniProt ID:	<u>P28472</u>
Synonyms:	DEE43; ECA5; EIEE43
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)