

Product datasheet for **HP200762**

GGCX Human qPCR Primer Pair (NM_000821)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2677
Forward Sequence:	CTCTACCTCCTGGAGCAGCTAT
Reverse Sequence:	GTCCCAGGAATAGCCATACAGC
ACCN:	NM_000821 , NM_000821.1 , NM_000821.2 , NM_000821.3 , NM_000821.4 , NM_000821.5 , NM_000821.6 , BC013979 , BC004422 , BM704752 , BM974022 , BQ000150 , NM_000821.7
UniProt ID:	P38435
Synonyms:	VKCFD1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)