

Product datasheet for **HP200778**

Metabotropic Glutamate Receptor 1 (GRM1) Human qPCR Primer Pair (NM_000838)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2911
Forward Sequence:	CTGGCATGAAGGAGTGCTGAAC
Reverse Sequence:	AGCAGCTCACTTCTCCTTCCG
ACCN:	<u>NM_000838</u> , <u>NM_000838.1</u> , <u>NM_000838.2</u> , <u>NM_000838.3</u> , <u>BC111844</u> , <u>BC136280</u> , <u>BC143777</u> , <u>BC143779</u>
UniProt ID:	<u>Q13255</u>
Synonyms:	GPRC1A; GRM1A; mGlu1; MGLUR1; MGLUR1A; SCAR13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)