

Product datasheet for **HP200789**

GSTM3 Human qPCR Primer Pair (NM_000849)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2947
Forward Sequence:	CGAAGCCAATGGCTGGATGTGA
Reverse Sequence:	GTTGTGCTTGCAGCGATGTAG
ACCN:	<u>BC030253</u> , <u>NM_000849</u> , <u>NM_000849.1</u> , <u>NM_000849.2</u> , <u>NM_000849.3</u> , <u>NM_000849.4</u> , <u>BC000088</u> , <u>BC008790</u> , <u>BE392471</u> , <u>BG285419</u> , <u>BT019945</u> , <u>NM_000849.5</u>
UniProt ID:	<u>P21266</u>
Synonyms:	GST5; GSTB; GSTM3-3; GTM3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)