

Product datasheet for **HP200794**

Glutathione S Transferase theta 2 (GSTT2) Human qPCR Primer Pair (NM_000854)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2953
Forward Sequence:	GAAACTGCCGACGCTCAAGGAT
Reverse Sequence:	GCAGGTCAGATGGATACCACTG
ACCN:	NM_000854 , NM_000854.1 , NM_000854.2 , NM_000854.3 , NM_000854.4 , BC002415.1 , BC146929 , BC146940 , BC146943 , BM804805 , BM976631 , BU584249
UniProt ID:	P0CG29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)