

## Product datasheet for **HP200803**

### HTR1B Human qPCR Primer Pair (NM\_000863)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3351
Forward Sequence:	CACCACGCTCTCCAATGCCTTT
Reverse Sequence:	CATCACCAGGATGGACACAAGC
ACCN:	<a href="#">NM_000863</a> , <a href="#">NM_000863.1</a> , <a href="#">NM_000863.2</a> , <a href="#">BC069065</a> , <a href="#">BC096206</a> , <a href="#">BC096207</a> , <a href="#">BC096208</a> , <a href="#">NM_000863.3</a>
UniProt ID:	<a href="#">P28222</a>
Synonyms:	5-HT-1B; 5-HT-1D-beta; 5-HT1B; 5-HT1DB; HTR1D2; HTR1DB; S12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)