

Product datasheet for **HP200818**

IL2 Receptor beta (IL2RB) Human qPCR Primer Pair (NM_000878)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3560
Forward Sequence:	GGTGAACCAAACCTGTGAGCT
Reverse Sequence:	GGTGACGATGTCAACTGTGGTC
ACCN:	NM_000878 , NM_000878.1 , NM_000878.2 , NM_000878.3 , NM_000878.4 , BC025691 , BC025691.1 , NM_000878.5
UniProt ID:	P14784
Synonyms:	CD122; IL15RB; IMD63; P70-75
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)