

## Product datasheet for **HP200821**

### IL12A Human qPCR Primer Pair (NM\_000882)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3592
Forward Sequence:	TGCCTTCACCACTCCCAAACCC
Reverse Sequence:	CAATCTCTCAGAAGTGCAAGGG
ACCN:	<u><a href="#">NM_000882</a></u> , <u><a href="#">NM_000882.1</a></u> , <u><a href="#">NM_000882.2</a></u> , <u><a href="#">NM_000882.3</a></u> , <u><a href="#">BC104982</a></u> , <u><a href="#">BC104982.1</a></u> , <u><a href="#">BC035571</a></u> , <u><a href="#">BC104984</a></u> , <u><a href="#">BG702253</a></u> , <u><a href="#">NM_000882.4</a></u>
UniProt ID:	<u><a href="#">P29459</a></u>
Synonyms:	CLMF; IL-12A; NFSK; NKSF1; P35
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)