

Product datasheet for **HP200823**

IMPDH2 Human qPCR Primer Pair (NM_000884)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3615
Forward Sequence:	AGTGGCTCCATCTGCATTACGC
Reverse Sequence:	GGATTCTCCATCAGCAATGACC
ACCN:	<u>NM_000884</u> , <u>NM_000884.1</u> , <u>BC006124</u> , <u>BC006124.1</u> , <u>BC009236</u> , <u>BC012840</u> , <u>BC015567</u> , <u>NM_000884.3</u>
UniProt ID:	<u>P12268</u>
Synonyms:	IMPD2; IMPDH-II
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)