

Product datasheet for **HP200838**

MGP Human qPCR Primer Pair (NM_000900)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4256
Forward Sequence:	CCTCAGCAGAGATGGAGAGCTA
Reverse Sequence:	ATGGCGTAGCGTTCGCAAAGTC
ACCN:	NM_000900 , NM_000900.1 , NM_000900.2 , NM_000900.3 , NM_000900.4 , BC093078 , BC005272 , BC070314 , BF699158 , BT006733 , NM_000900.5
UniProt ID:	P08493
Synonyms:	GIG36; MGLAP; NTI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)