

Product datasheet for **HP200846**

Natriuretic Peptide Receptor C (NPR3) Human qPCR Primer Pair (NM_000908)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4883
Forward Sequence:	CAGTGGAGACTACGCCTTCTTC
Reverse Sequence:	TGACTGTCTGGAGGGACGAGTA
ACCN:	<u>NM_000908</u> , <u>NM_000908.1</u> , <u>NM_000908.2</u> , <u>BC022250</u> , <u>BC131540</u> , <u>BF224187</u> , <u>BQ014050</u> , <u>NM_000908.4</u>
UniProt ID:	<u>P17342</u>
Synonyms:	ANP-C; ANPR-C; ANPRC; C5orf23; GUCY2B; NPR-C; NPRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)