

## Product datasheet for **HP200859**

### PDE3A Human qPCR Primer Pair (NM\_000921)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5139
Forward Sequence:	AGAGCCTCTGAGGAAAGCATCG
Reverse Sequence:	GCTGCTCCATAATTGAGTCCAGG
ACCN:	<u><a href="#">NM_000921</a></u> , <u><a href="#">NM_000921.1</a></u> , <u><a href="#">NM_000921.2</a></u> , <u><a href="#">NM_000921.3</a></u> , <u><a href="#">NM_000921.4</a></u> , <u><a href="#">BC117369</a></u> , <u><a href="#">BC117371</a></u> , <u><a href="#">BQ024709</a></u> , <u><a href="#">BQ187047</a></u> , <u><a href="#">BX646774</a></u> , <u><a href="#">NM_000921.5</a></u>
UniProt ID:	<u><a href="#">Q14432</a></u>
Synonyms:	CGI-PDE; CGI-PDE-A; CGI-PDE A; HTNB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)