

Product datasheet for **HP200870**

Phospholipase C beta 4 (PLCB4) Human qPCR Primer Pair (NM_000933)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5332
Forward Sequence:	CTGACCAAGAGGAGGAAGCTCA
Reverse Sequence:	CCTGCTCATCTTCTACAGTGACC
ACCN:	NM_000933 , NM_000933.1 , NM_000933.2 , NM_000933.3 , BC117458 , BC143868 , BG220268 , BI439140 , BX648680 , NM_000933.4
UniProt ID:	Q15147
Synonyms:	ARCND2; PI-PLC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)