

Product datasheet for **HP200873**

Pancreatic Lipase (PNLIP) Human qPCR Primer Pair (NM_000936)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5406
Forward Sequence:	CACACAAGCCTCGCAGAACATC
Reverse Sequence:	TGGAAGGTGAGTAACCGAACGC
ACCN:	<u>BC014309</u> , <u>NM_000936</u> , <u>NM_000936.1</u> , <u>NM_000936.2</u> , <u>NM_000936.3</u> , <u>BC014309.2</u> , <u>BC007025</u>
UniProt ID:	<u>P16233</u>
Synonyms:	PL; PNLIPD; PTL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)