

## Product datasheet for **HP200874**

### **POLR2A Human qPCR Primer Pair (NM\_000937)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	5430
<b>Forward Sequence:</b>	GAGAGCGTTGAGTCCAGAACC
<b>Reverse Sequence:</b>	TGGATGTGTGCGTTGCTCAGCA
<b>ACCN:</b>	<u><a href="#">NM_000937</a></u> , <u><a href="#">NM_000937.1</a></u> , <u><a href="#">NM_000937.2</a></u> , <u><a href="#">NM_000937.3</a></u> , <u><a href="#">NM_000937.4</a></u> , <u><a href="#">BC067295</a></u> , <u><a href="#">BC107140</a></u> , <u><a href="#">BC137231</a></u> , <u><a href="#">BQ353795</a></u> , <u><a href="#">BQ891728</a></u> , <u><a href="#">BU730469</a></u> , <u><a href="#">BU739765</a></u> , <u><a href="#">NM_000937.5</a></u>
<b>UniProt ID:</b>	<u><a href="#">P24928</a></u>
<b>Synonyms:</b>	hRPB220; hsRPB1; NEDHIB; POLR2; POLRA; RPB1; RPBh1; RpIILS; RPO2; RPOL2
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)