

## Product datasheet for **HP200886**

### Prolactin Receptor (PRLR) Human qPCR Primer Pair (NM\_000949)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5618
Forward Sequence:	CATGGTGACCTGCATCTTTCCG
Reverse Sequence:	GTGGGAGGAAAGTCTTGGCATC
ACCN:	<a href="#">NM_000949</a> , <a href="#">NM_000949.1</a> , <a href="#">NM_000949.2</a> , <a href="#">NM_000949.3</a> , <a href="#">NM_000949.4</a> , <a href="#">NM_000949.5</a> , <a href="#">NM_000949.6</a> , <a href="#">BC059392</a> , <a href="#">BC059392.1</a> , <a href="#">BP357907</a> , <a href="#">BQ023302</a> , <a href="#">BX538149</a> , <a href="#">NM_000949.7</a>
UniProt ID:	<a href="#">P16471</a>
Synonyms:	HPRL; hPRLrI; MFAB; RI-PRLR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)