

Product datasheet for **HP200887**

PRRG1 Human qPCR Primer Pair (NM_000950)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5638
Forward Sequence:	AGACAGGGCAACATTGAGCGTG
Reverse Sequence:	GCGCTTTTGTGTAGGTGCTCCA
ACCN:	<u>NM_000950</u> , <u>NM_000950.1</u> , <u>NM_000950.2</u> , <u>BC060833</u> , <u>BC060833.1</u> , <u>BC012608</u> , <u>BC030786</u> , <u>BC041591</u> , <u>BX097465</u> , <u>NM_000950.3</u>
UniProt ID:	<u>O14668</u>
Synonyms:	PRGP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)