

Product datasheet for **HP200895**

EP4 (PTGER4) Human qPCR Primer Pair (NM_000958)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5734
Forward Sequence:	TACTCATTGCCACCTCCCTGGT
Reverse Sequence:	GACTTCTCGCTCCAAACTTGGC
ACCN:	<u>NM_000958</u> , <u>NM_000958.1</u> , <u>NM_000958.2</u> , <u>BC101534</u> , <u>BC101534.1</u> , <u>BC027934</u> , <u>BC113523</u> , <u>BM312677</u> , <u>BQ029720</u> , <u>NM_000958.3</u>
UniProt ID:	<u>P35408</u>
Synonyms:	EP4; EP4R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)