

Product datasheet for **HP200899**

PTGS1 Human qPCR Primer Pair (NM_000962)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5742
Forward Sequence:	GATGAGCAGCTTTTCCAGACGAC
Reverse Sequence:	AACTGGACACCGAACAGCAGCT
ACCN:	NM_000962 , NM_000962.1 , NM_000962.2 , NM_000962.3 , BC029840 , BC029840.1
UniProt ID:	P23219
Synonyms:	COX1; COX3; PCOX1; PES-1; PGG/HS; PGHS-1; PGHS1; PHS1; PTGHS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)