

Product datasheet for **HP200900**

COX2 (PTGS2) Human qPCR Primer Pair (NM_000963)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5743
Forward Sequence:	CGGTGAACTCTGGCTAGACAG
Reverse Sequence:	GCAAACCGTAGATGCTCAGGGA
ACCN:	<u>NM_000963</u> , <u>NM_000963.1</u> , <u>NM_000963.2</u> , <u>NM_000963.3</u> , <u>BC013734</u> , <u>BC013734.1</u> , <u>BM978796</u> , <u>NM_000963.4</u>
UniProt ID:	<u>P35354</u>
Synonyms:	COX-2; COX2; GRIPGHS; hCox-2; PGG/HS; PGHS-2; PHS-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)