

Product datasheet for **HP200914**

RPL23 Human qPCR Primer Pair (NM_000978)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9349
Forward Sequence:	ATCAAGGGACGGCTGAACAGAC
Reverse Sequence:	GTCGAATGACCACTGCTGGATG
ACCN:	BC034378 , NM_000978 , NM_000978.1 , NM_000978.2 , NM_000978.3 , BC104651 , BC104651.1 , BC003518 , BC010114 , BC062716 , BC106061 , NM_000978.4
UniProt ID:	P62829
Synonyms:	L23; rpL17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)