

## Product datasheet for **HP200928**

### RPL29 Human qPCR Primer Pair (NM\_000992)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 6159   |
| Forward Sequence: | CCCGATCACAAAGATACGAATCTC   |
| Reverse Sequence: | TGCACTCATGGCCTTGCCATTG   |
| ACCN:             | <u><a href="#">NM_000992</a></u> , <u><a href="#">NM_000992.1</a></u> , <u><a href="#">NM_000992.2</a></u> , <u><a href="#">BC008926</a></u> , <u><a href="#">BC018943</a></u> , <u><a href="#">BC070190</a></u> , <u><a href="#">BC070481</a></u> ,<br><u><a href="#">BC071663</a></u> , <u><a href="#">BC071909</a></u> , <u><a href="#">NM_000992.3</a></u> |
| UniProt ID:       | <u><a href="#">P47914</a></u>  |
| Synonyms:         | HIP; HUMRPL29; L29; RPL29P10; RPL29_3_370  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                   |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)