

Product datasheet for **HP200937**

RPL36AL Human qPCR Primer Pair (NM_001001)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6166
Forward Sequence:	ATCGGAAGCAGAGTGGCTATGG
Reverse Sequence:	CAGCATCCTCTTGATCTGCAG
ACCN:	NM_001001 , NM_001001.1 , NM_001001.2 , NM_001001.3 , NM_001001.4 , BC107716 , BC107716.1 , BC000741 , BC003145 , BC070207 , BM973694 , NM_001001.5
UniProt ID:	Q969Q0
Synonyms:	RPL36A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)