

Product datasheet for **HP200946**

REEP3 Human qPCR Primer Pair (NM_001001330)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221035
Forward Sequence:	CAGCTACTGCTGCTGTTACTGC
Reverse Sequence:	GGTCTTTGTCCCACAGGCTCAT
ACCN:	<u>BC018658</u> , <u>NM_001001330</u> , <u>NM_001001330.1</u> , <u>NM_001001330.2</u> , <u>BC018658.1</u> , <u>BC068557</u> , <u>BC010040</u> , <u>BC038724</u> , <u>BC057832</u> , <u>NM_001001330.3</u>
UniProt ID:	<u>Q6NUK4</u>
Synonyms:	C10orf74; Yip2b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)